

## **MediSelect II Oxygen Piston Regulators**



## Features:

- Regulator with Flow Selector
- Rotating pressure gauge which allows convenient reading
- 360° swivelling outlet it enables better orientation of the nasal cannula or oxygen mask towards the patient (preventing from twisting)
- Innovative self centering flow setting device with continuous flow between settings. In the unlikely event of indent mechanism failure, the patient will still be supplied by medicinal gas
- Lateral and frontal reading of flow settings
- Higher number of flow disc holes increases treatment options. Extra flow setting of 25I/min on the traditional 15I/min variant, allows use in resuscitation

## Part Numbers:

551189	Oxygen Piston Regulator 0-25L/min, Barbed Nipple, two S.I.S Sleeves (1, 2, 3, 4, 5, 6, 7, 9, 12, 15 & 25L/min)
551191	Oxygen Piston Regulator 0-15L/min, Barbed Nipple (1, 2, 3, 4, 5, 6, 7, 9, 12 &15L/min)
551197	Oxygen Piston Regulator 0-2L/min, Barbed Nipple (Low Flow) (0.1, 0.2, 0.3, 0.4, 0.5, 0.6, 0.7, 0.8, 1.0, 1.5, 2L/min)



w w w . parkerhealth.com.au w w w . parkerhealth.com . m y w w w . parkerhealth.co.nz



Scan QR code to download pdf

PD 941